


Possibilities for Studying Sulfur Speciation in Biological Systems Using MEX



XAS & MEX Workshop, August 26th 2021
Dr Mark Hackett – Curtin University



Curtin University

BACKGROUND

- ❖ BSc (Hons) Analytical Chemistry (Curtin University, 2006)
- ❖ PhD Chemistry (The University of Sydney, 2007 - 2010)
 - *Spectroscopic investigation of cerebral malaria*
 - *ischemia and blood-brain-barrier*
- ❖ Research Fellowships in Canada (University of Saskatchewan 2011 - 2015)
 - *Spectroscopic studies of neurodegeneration after stroke*
 - *ischemic stroke, intra-cerebral haemorrhage*
- ❖ Senior Lecturer in Analytical Chemistry at Curtin University (2016)
 - *Developing and applying spectroscopic approaches to study chemical mechanisms of memory loss during ageing and dementia*



Analytical Chemistry
Applied Spectroscopy

Neuroscience
Neurochemistry



CURTIN HEALTH
INNOVATION
RESEARCH INSTITUTE



Curtin University



Australian Government

Australian Research Council

CHEMISTRY @ CURTIN





TALK OUTLINE



- Application of XAS to image Zn speciation in brain tissue
- Development of chemically specific XANES sulfur imaging
- Pitfalls with chemically specific sulfur imaging
- Applications of chemically specific sulfur imaging to study brain tissue

SAMPLE PREPARATION – RAPID FLASH FREEZING


- Tissues samples are rapidly dissected and immediately flash frozen in liquid-nitrogen cooled iso-pentane
 - *(no perfusion, no fixation, no sucrose cryoprotection, no plastic or paraffin embedding)*
- Samples stored in -80°C until ~ 1 week before analysis
- Samples sectioned on a *cryo-microtome at -16°C* for brain (liver -18°C), *$10\text{-}\mu\text{m-thick}$* . This is warmer than most SOPs, which are written for sucrose-cryoprotected, formalin-fixed tissue
- Tissue sections “melted” onto substrate, and either *transferred back to -80°C freezer, or air-dried*
- Analysis typically within 1 week of sectioning
- If your samples are cracking, or tearing when you section, they are too cold. It takes about 30 – 40 min for tissue to warm from -80°C to -16°C in cryostat
- Substrates
 - 500 nm, 1 μm , or 2 μm thick Si_3N_4
 - *XFM, or multi-modal XFM/IRM*
 - 500 μm or 1 mm thick CaF_2 , BaF_2
 - *For IRM*



RAPID FLASH FREEZING, HISTOLOGY, & ICE CRYSTALS




- Histology and immuno-histochemistry can easily be performed on flash frozen tissues
- Either on the same tissue analysed at synchrotron, or serial sections
- Post-sectioning immersion fixation (typically 10 % PFA)

Freezing on dry ice




Neuroimage. 2011, 55(1): 32–38.

Liq. N₂ cooled isopentane







Liq. N₂ cooled isopentane H&E




DEVELOPING XANES Imaging – Zn in Brain Tissue

XANES spectroscopy is highly sensitive to the coordination geometry, number of ligands, and ligand type




IS THE ZINC K-EDGE SENSITIVE TO COORDINATION CHEMISTRY? (YES)

Modelling Zn Bonding in Biological Environments: Zn^{2+} in presence of excess biological ligands





ZN SPECIATION IN BRAIN TISSUE





Ashley Hollings PhD

Hollings et. al., *Metallomics*, 12, 2143-2144, 2020.



CHEMICALLY SPECIFIC XANES ZN IMAGING (XFM)



~1 ms dwell
~ 2 x 2 mm area
~2 – 5 μm pixel
~ 8 hours
~1 hours (ROI)




A CAUTIONARY TALE: SAMPLE PREPARATION & ANALYSIS CONDITIONS ARE CRITICAL!!



Zn


Max Min

Scale = 500 μ m




Lee, Joo-Yong, et. al., *Brain research* 1418, 12-22, **2011**





Hollings et. al., *JAAS*, 35, 2498-2508. **2020**




A CAUTIONARY TALE: SAMPLE PREPARATION & ANALYSIS CONDITIONS ARE CRITICAL!!



**Experimental considerations for the XANES library:
buffer ions, and “glassing agents”**





MEDIUM ENERGY X-RAY ABSORPTION SPECTROSCOPY (MEX) @ THE SYNCHROTRON



MEDIUM ENERGY X-RAY ABSORPTION SPECTROSCOPY

The Sulfur K-edge Is Very Sensitive To Sulfur Speciation
Fitting Sulfur K-edge XANES of Brain Tissue



BE WARY OF SAMPLE THICKNESS, AND CONCENTRATION

- Self-Absorption distorts XANES spectra
- Fitting errors
- Self-absorption is often not negligible in medium-energy range
- We aim for thin samples (< 20 μm)
- Aqueous standard solutions
- CXRO Website
- https://henke.lbl.gov/optical_constants/

Research article

Received: 11 April 2012

Accepted: 31 May 2012


Published online in Wiley Online Library: 29 June 2012

(wileyonlinelibrary.com) DOI 10.1002/xrs.2407


X-RAY
Spectrometry

Sample thickness considerations for quantitative X-ray fluorescence analysis of the soft and skeletal tissues of the human body – theoretical evaluation and experimental validation

Magdalena Szczerbowska-Boruchowska*



Pickering et. al., FEBS Letters 2008






Pickering et. al., Biochemistry 2001


MEDIUM ENERGY X-RAY ABSORPTION SPECTROSCOPY TIME CONSTRAINTS

- Typically less flux at MEX beamlines
 - (relative to hard X-ray microprobes)
- Dwell time (50 – 500 ms)
- Full XANES per pixel, generally not feasible (time)

$$I_e = A_{\text{sulfonic acid}} \times I_{N(\text{sulfonic acid})} + A_{\text{sulfate ester}} \times I_{N(\text{sulfate ester})}$$







PIT FALLS WITH CHEMICALLY SPECIFIC XFM





Energy Calibration & Stability is Critical!

X-RAY PHOTO-DAMAGE



SAMPLE PREPARATION – DON'T USE CHEMICAL FIXATION



Chemical fixation = sulfonic acid leaching & thiol oxidation
Paraffin or plastic embedding = massive sulfur oxidation



Hackett et al. (2012) ACS Chemical Neuroscience


Hackett et al. (2016) ACS Biochemistry

EFFECTS OF SAMPLE PREPARATION






EFFECTS OF POST-MORTEM INTERVAL

The time it takes to sacrifice an animal and remove its brain (30 seconds – 3 minutes) induces ischemia and oxidises thiols to disulfides



MEDIUM ENERGY X-RAY ABSORPTION SPECTROSCOPY

BRAIN TAURINE LEVELS AFTER STROKE




MEX: IS THERE AN ISCHEMIC PENUMBRA IN THE MURINE PHOTOTHROMBOTIC STROKE MODEL?

H&E



GFAP




Scale = 500 μ m







Saskatchewan Cerebrovascular Unit

Rose Bengal



Is there an ischemic penumbra in the photothrombotic stroke model?

MEX: IS THERE AN ISCHEMIC PENUMBRA IN THE MURINE PHOTOTHROMBOTIC STROKE MODEL?



ACKNOWLEDGEMENTS



Curtin University

Curtin

Meg Willans
Ashley Hollings
Karina Khambatta
David Hartnell
Rhianon Boseley
Gae Ellison

USASK

Mike Kelly
Phyllis Paterson
Graham George
Ingrid Pickering
Jake Pushie

Curtin – CHIRI

John Mamo
Ryu Takechi
Virginie Lam
Ross Graham
Clinton Kidman

UWA

Peter Arthur
Jessica Terrill
Nicole Smith

Monash

Ramesh Rajan

SSRL

Sam Webb

ANSTO

David Paterson
Martin de Jonge
Daryl Howard
Michael Jones
Cameron Kewish
Mark Tobin
Keith Bambery
Pimm Vongsvivut
Annaleise Klein
Danielle Martin
Reiner Siegele
David Cohen



CURTIN HEALTH
INNOVATION
RESEARCH INSTITUTE







Australian Government

Australian Research Council



INTRODUCTION: THE TECHNIQUES I USE.....

XRF (XFM)



WHAT IS THE ROLE OF SPECTROSCOPIC IMAGING?

